



MEMO

Valued Customer

To

Consultation & Reference Laboratory

From

Washed Red Blood Cell Product Codes – ACP 215

Subject

6/20/2024

Date

Dear Valued Partner,

In approximately one month, our Consultation and Reference laboratory will be manufacturing washed red blood cells (RBCs) on the ACP 215 automated cell processor. RBCs processed on the ACP 215 are prepared in a functionally closed system. The product expiration of washed RBCs on the ACP 215 is 24 hours from the date and time the unit is prepared.

The product codes are as follows:

- E5638 Washed RED BLOOD CELLS|None/XX/refg|ResLeu:<5E6
- E5639 Washed RED BLOOD CELLS|None/XX/refg|Irradiated|ResLeu:<5E6

Please let me know if you have any questions.

Best Regards,

A handwritten signature in cursive script that reads "Shantelle Myers".

Shantelle Myers, MLS(ASCP)SBB^{CM}
Manager, Consultation & Reference Laboratory
Gulf Coast Regional Blood Center
(713) 791-6278
smyers@giveblood.org

Commit for Life.[®]

1400 La Concha Lane, Houston, Texas 77054 • 713-790-1200 • 1-888-482-5663 • www.giveblood.org

A Nonprofit Community Blood Center and Member of America's Blood Centers, aaBB and the Texas Medical Center.